


Characterization of Novel PDMS-Based Drug Encapsulation Strategies for Controlled Release Applications Using Choline-Geranic Acid Ionic Liquid

Aasim Khan¹, Rachel Reyes De Los Santos², Mallikarjun Gosangi^{3,4}, Kaushal Rege^{3,4,5}

¹Biomedical Engineering Graduate Program, School of Biological and Health Systems Engineering, Arizona State University, Tempe, Arizona, USA
²Biochemistry Undergraduate Program, School of Molecular Sciences, Arizona State University, Tempe, Arizona, USA
³Center for Biomaterials Innovation and Translation, Biodesign Institute, Arizona State University, Tempe, Arizona, USA
⁴Chemical Engineering, School for Engineering of Matter, Transport, and Energy, Arizona State University, Tempe, Arizona, USA
⁵Biological Design Graduate Program, School for Engineering of Matter, Transport, and Energy, Arizona State University, Tempe, Arizona, USA

Introduction


Polydimethylsiloxane (PDMS) is an elastomer commonly used in biomedical implants due to its favorable mechanical properties and biocompatibility.



- Combining PDMS with a drug delivery platform enhances its biomedical utility by providing both structural function and therapeutic benefit
- We developed a novel formulation combining PDMS with ionic liquids that enables high loading of the model drug Rhodamine B and supports sustained, high-level release over time.

Ionic Liquids


¹H NMR Spectrum of Ionic liquid



¹H-NMR spectroscopy of CAGE ionic liquids made up of four different molar ratios of cation, Choline, and anion, Geranic Acid, molecules.


Methods

Fabrication of RhB-loaded PDMS devices




PDMS films were prepared by mixing Sylgard 184 base with curing agent at a 10:1 ratio.

Drug delivery as films




Drug delivery as implants



Results

SEM Imaging



Loading Quantities

| Average weight and RhB Loading into films | | | |
|---|---------------------|-----------------|-------------|
| Thin film (1 mm x 10 mm) | Wt of the film (mg) | RhB loaded (µg) | RhB (µg/mg) |
| PDMS+DCM+RhB | 98.6 | 90.4 | 0.9 |
| PDMS+H2O+RhB | 81.6 | 69.6 | 0.9 |
| PDMS+IL-2+RhB | 89.1 | 225.3 | 2.5 |

| Average weight and RhB Loading into Implants | | | |
|--|---------------------|-----------------|-------------|
| Thick film (7 mm x 10 mm) | Wt of the film (mg) | RhB loaded (µg) | RhB (µg/mg) |
| PDMS+DCM+RhB | 538.4 | 335.6 | 0.6 |
| PDMS+H2O+RhB | 550.4 | 335.7 | 0.6 |
| PDMS+IL-2+RhB | 615.8 | 1093.0 | 1.8 |


Korsmeyer–Peppas Model

| Film Type | Release Constant (k) | Release Exponent (n) |
|------------|----------------------|----------------------|
| DCM | 0.07230 | 0.36730 |
| H2O | 0.08070 | 0.47322 |
| CAGE (1:1) | 0.33747 | 0.36234 |
| CAGE (1:2) | 0.37259 | 0.29446 |
| CAGE (1:3) | 0.40952 | 0.28033 |
| CAGE (1:4) | 0.43352 | 0.26854 |

| Film Type | Release Constant (k) | Release Exponent (n) |
|------------|----------------------|----------------------|
| DCM | 0.01769 | 0.51708 |
| H2O | 0.01708 | 0.58465 |
| CAGE (1:2) | 0.19946 | 0.40080 |


Mechanical Testing

Young Modulus



- CAGE incorporation significantly reduced stiffness (~0.3–0.5 MPa), compared to 1:5 PDMS (~2.4 MPa).
- Lower modulus supports flexibility for coating applications.


Tensile Strength




- CAGE films had reduced tensile strength (< 2 MPa) relative to 1:5 PDMS (> 4 MPa).
- Mechanical properties remain suitable for non-load-bearing uses.

RhB Release from devices


Release from Film




- CAGE films released up to ~150 µg RhB, ~3× more than DCM or water.
- IL-1 to IL-3 reached >90% release; DCM and water stayed below 50%.
- k increased with CAGE ratio (up to 0.43), showing enhanced release rates.
- n decreased with CAGE ratio (down to 0.27), indicating Fickian diffusion.



Release from Implant



- CAGE films released up to ~1500 µg RhB; DCM and water stayed <350 µg.
- IL-2 released ~75% of loaded RhB; others remained <25%.
- k for IL-2 was 0.20 — nearly 10× higher than DCM or water films.
- n for IL-2 was 0.40, suggesting predominantly Fickian transport.



Conclusions / Future Directions

- CAGE-enhanced PDMS films enabled tunable drug loading and sustained release, offering a promising platform for localized delivery in implant coatings and reconstructions.
- Future directions include applying this strategy toward mitigating foreign body response (FBR) using small molecule therapeutics in long-term implants.

REFERENCES

- Marmo, A.C., Grunlan, M.A., Biomedical Silicones: Leveraging Additive Strategies to Propel Modern Utility. *ACS Macro Lett.* <https://doi.org/10.1021/acsmacrolett.2c00701>. (2023)
- Miranda et al. Properties and Applications of PDMS for Biomedical Engineering: A Review. *Journal of Controlled Release* <https://pubs.rsc.org/en/article/PMC8788510>. (2021)
- Padmanabhan et al. Allometrically Scaling tissue forces drive pathological foreign-body responses to implants via Rac2-activated myeloid cell. *Nat. Biomed. Eng.* <https://www.nature.com/articles/s41551-023-01091-5> (2023)

Acknowledgements

I sincerely thank Dr. Kaushal Rege, principal investigator of the Center for Biomaterials Innovation and Translation at Arizona State University, for his guidance and support throughout this project. I am especially grateful to Mallikarjun Gosangi for his mentorship, technical insights, and constant encouragement during every stage of this research. I also appreciate the resources and collaborative environment provided by the CBIT Lab, which were essential to the successful completion of this work.